




# Recombinant Mumps virus Hemagglutinin-neuraminidase (HN), partial

<b>Product Code</b>	CSB-EP319979MJJ2
<b>Abbreviation</b>	Recombinant Mumps virus Hemagglutinin-neuraminidase protein, partial
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P11235
<b>Storage Buffer</b>	Tris-based buffer,50% glycerol
<b>Product Type</b>	Recombinant Proteins
<b>Immunogen Species</b>	Mumps virus (strain Miyahara vaccine) (MuV)
<b>Purity</b>	≥ 85% as determined by SDS-PAGE.
<b>Sequence</b>	YATHDFSIGHPLNMPFSIPTATSPNGCTRIPSFSLGKTHWCYTHNVINANCKDH TSSNQYISMGILVQTASGYPMFKTLKIQYLSDGLNRKSCSIATVPDGCAMYCYV STQLETDDYAGSSPPTQKLTLFFYNQVTVTERTISPTGLEGNWATLVPGVGSGIY FENKLIFPAYGGVLPNSSLGVKSAREFFRPVNPYNPCSGPQQDLQDR
<b>Research Area</b>	Epigenetics and Nuclear Signaling
<b>Source</b>	E.coli
<b>Target Names</b>	HN
<b>Protein Names</b>	Recommended name: Hemagglutinin-neuraminidase EC= 3.2.1.18
<b>Expression Region</b>	152-360aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal 6xHis-tagged
<b>Mol. Weight</b>	27.0 kDa
<b>Protein Length</b>	Partial
<b>Image</b>	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at  $-20^{\circ}\text{C}/-80^{\circ}\text{C}$ . The shelf life of lyophilized form is 12 months at  $-20^{\circ}\text{C}/-80^{\circ}\text{C}$ .