




Recombinant Mouse Desmoglein-3 (Dsg3), partial

Product Code	CSB-EP007205MOb0
Abbreviation	Recombinant Mouse Dsg3 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O35902
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Mus musculus (Mouse)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	EWVKFAKPCREREDNSRRNPIAKITSDFQKNQKITYRISGVGIDQPPFGIFVVD PNNGDINITAIVDREETPSFLITCRALNALGQDVERPLILTVKILDVNDNPPIFSQT IFKGEIEENSASNSLVMILNATDADEPNHMNSKIAFKIVSQEPAGMSMFLISRNT GEVRTLTSSLDREQISSYHLVVSGADNDGTGLSTQCECSIKIKDVNDNFPVLR SQYSARIEENTLNAELLRFQVTDWDEEYTDNWLAVYFFFTSGNEGNWFEIETDP RTNEGILKVVKALDYEQVQSMQFSIAVRNKAEFHQSVISQYRVQSTPVTIQVID VREGISFRPPSKTFTVQRGVSTNKLVG YILGTYQATDEDTGKAASSVRYVLGR NDGGLLVIDSKTAQIKFVKNIDRDSTFIVNKTISAEVLAIDENTGKTSTGTIYVEVP SFNENCPSVVLEKKDICTSSPSVTLVSVRTLDRGKYTGPLYTVSLEEQLKLPVM WTITTLNATSALLQAQQQVSPGVYNVPVIVKDNQDGLCDTPESLTLTVQCDD RSMCRAPIPSREPNTYGESSWRLGP
Research Area	Others
Source	E.coli
Target Names	Dsg3
Expression Region	50-617aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged
Mol. Weight	66.9 kDa
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.