




Recombinant Arachis hypogaea Allergen Ara h 1, clone P17, partial

Product Code	CSB-EP331779ANEa0
Abbreviation	Recombinant Arachis hypogaea Allergen Ara h 1, clone P17 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P43237
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Arachis hypogaea (Peanut)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	KSPYRK TENPCAQRCLQSCQQEPDDLKQKACESRCKLEYDPRCVYDTGAT NQRHPPGERTRGRQPGDYDDRRQRREEGGRWGPAEPREREREEDWRQ PREDWRRPSHQQPRKIRPEGREGEQEWGTPGSEVREETSRRNPFYFPSRRF STRYGNQNGRIRVLQRFQDQRSKQFQNLQNHRIQIEARPN
Research Area	Others
Source	E.coli
Target Names	N/A
Expression Region	26-216aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	29.0 kDa
Protein Length	Partial

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.