




Recombinant Mycobacterium paratuberculosis UPF0353 protein MAP_3434 (MAP_3434), partial

Product Code	CSB-EP748202MLA1
Abbreviation	Recombinant Mycobacterium paratuberculosis UPF0353 protein MAP_3434, partial protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q73UD4
Storage Buffer	Tris-based buffer, 50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Mycobacterium paratuberculosis (strain ATCC BAA-968 / K-10)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	NRAVIMLVIDMSQSMRATDVEPNRLKAAEQAASQFASQLTPGINLGLVGFAGT PYLLVPPTPQHQTIDALKKLDFASTATGEAIFTALHAISATAVAGGDTPPPAR IVLLSDGGENKPSNPSPDHPDGVYTAARLAKDEGVPISTITFGTKGGEIEMDGQK VAVPVSTDQMKMVAKLSSGGQSYTATNLGELQKSYNAIENE
Research Area	Others
Source	E.coli
Target Names	MAP_3434
Protein Names	Recommended name: UPF0353 protein MAP_3434
Expression Region	91-292aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged
Mol. Weight	27.1 kDa
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.