




# Recombinant Rabies virus Glycoprotein (G), partial

<b>Product Code</b>	CSB-EP365970RIDa0
<b>Abbreviation</b>	Recombinant Rabies virus Glycoprotein, partial
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P03524
<b>Storage Buffer</b>	Tris-based buffer,50% glycerol
<b>Product Type</b>	Recombinant Proteins
<b>Immunogen Species</b>	Rabies virus (strain ERA)(RABV)
<b>Purity</b>	≥ 90% as determined by SDS-PAGE.
<b>Sequence</b>	KFPIYTILDKLGPWSPIDIHHLSCPNNLVEDEGCTNLSGFSYMEKLVGYILAIAIK MNGFTCTGVVTEAETYTNFVGIVTTTTFKRKHFRPTPDACRAAYNWKMAGDP RYEESLHNPYPDYRWLRTVKTTKESLVIISPSVADLDPYDRSLHSRVFSPGKCS GVAVSSTYCSTNHDIWMPENPRLGMSCDIFTNSRGKRASKGSETCGFVDE RGLYKSLKGACKLKLGVGLRLMDGTWVAMQTSNETKWCPPDQLVNLHDF RSDEIEHLVVEELVRKREECLDALESIMTTKSVSFRRLSHLRKLVPGFVKAYTIF NKTLM EADAHYKSVRTWNEILPSKGCLR VGG RCHPHVNGVFFNGIILGPDGNV LIPEMQSSLLQQHMELLESSVIPLVHPLADPSTVFKDGD EAE DFVEVHLPDVHN QVSGVDLGLPNWGKY
<b>Research Area</b>	Signal Transduction
<b>Source</b>	E.coli
<b>Target Names</b>	G
<b>Expression Region</b>	20-459aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal 6xHis-tagged
<b>Mol. Weight</b>	55.0 kDa
<b>Protein Length</b>	Partial
<b>Image</b>	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.