




Recombinant Staphylococcus aureus Enterotoxin type A (entA), partial

Product Code	CSB-EP362952FKZ1
Abbreviation	Recombinant Staphylococcus aureus entA protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P0A0L2
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Staphylococcus aureus
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	SEKSEEINEKDLRKKSELQGTALGNLQIYYNEKAKTENKESHQFLQHTILF KGFFTDHSWYNDLLVDFDSKDIVDKYKGGKVDLYGAYGYQCAGGTPNKTAC MYGGVTLHDNNRLTEEKVPINLWLDGKQNTVPLETVKTNKKNVTQELDLQA RRYLQEKYNLYNSDVFDGKVRGLIVFHTSTEPSVNYDLFGAQQQYSNTLLRI YRDNKTINSENMHIDIY
Research Area	Cell Biology
Source	E.coli
Target Names	entA
Expression Region	25-253aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	30.5
Protein Length	Partial

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.