




# Recombinant Mouse Angiopoietin-like protein 8 (Angptl8)

<b>Product Code</b>	CSB-BP844436MO
<b>Abbreviation</b>	Recombinant Mouse Angptl8 protein
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8R1L8
<b>Storage Buffer</b>	Tris-based buffer,50% glycerol
<b>Product Type</b>	Recombinant Proteins
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	≥ 85% as determined by SDS-PAGE.
<b>Sequence</b>	VRPAPVAPLGGPEPAQYEELTLLFHGALQLGQALNGVYRATEARL TEAGHSLG LYDRALEFLGTEVRQGQDATQELRTSLSEIQVEEDALHLRAEATARSLGEVAR AQQALRDTVRRLQVQLRGAWLGQAHQEFETLKARADKQSHLLWALTGHVQR QQREMAEQQQWLRQIQQLHTAALPA
<b>Research Area</b>	Cell Biology
<b>Source</b>	Baculovirus
<b>Target Names</b>	Angptl8
<b>Protein Names</b>	Recommended name: Hepatocellular carcinoma-associated protein TD26 homolog
<b>Expression Region</b>	16-198aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal 10xHis-tagged
<b>Mol. Weight</b>	23.0 kDa
<b>Protein Length</b>	Full Length of Mature Protein
<b>Image</b>	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.