




Recombinant Mouse Receptor-interacting serine/threonine-protein kinase 1 (Ripk1)

Product Code	CSB-BP720181MO
Abbreviation	Recombinant Mouse Ripk1 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q60855
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Mus musculus(Mouse)
Purity	Greater than $\geq 85\%$ as determined by SDS-PAGE.
Sequence	MQPDMSLDNIKMASSDLLEKTDLDSSGGFGKVS LCYHRSHGFVILKKVYTGPNR AEYNEVLLEEGKMMHRLRHSRVVKLLGIIIEEGNYSLVMEYMEKGNLMHV LKT QIDVPLSLKGR IIVEAIEGMCYLHDKGV IHKDLKPENILVDRDFHIKIADLG VASFK TWSKLTKEKDNKQKEV SSTTKKNGGTLYYMAPEHLNDINAKPTEKSDVYSF GIVLWAIFAKKEPYENVICTEQFVICIKSGNRPNVEEILEYCPREIISLMERCWQA IPEDRPTFLGIEEEFRPFYLSHFEEYVEEDVASLKKEYPDQSPVLQRMFSLQHD CVPLPPSRSNSEQPGSLHSSQGLQMGPVEESWFSSSPEYPQDENDRSVQAK LQEEASYHAFGIFA EKQTKPQPRQNEAYNREEERKRRVSHDPFAQQRARENI KSAGARGHSDPSTTSRGI AVQQLSWPATQTVWNNGLYNQHGF GTTGTGVVY PPNLSQMYSTYKTPVPETNIPGSTPTMPYFSGPVADDLIKYTIFNSSGIQIGNHN YMDVGLNSQPPNNTCKEESTSRHQ AIFDNTTSLTDEHLNPIRENLGRQWKNC ARKLGFTESQIDEIDHDYERDGLKEKVYQMLQKWL MREGTKGATVGKLAQAL HQCCRIDLLNHLIRASQS
Research Area	Cancer
Source	Baculovirus
Target Names	Ripk1
Protein Names	Recommended name: Receptor-interacting serine/threonine-protein kinase 1 EC= 2.7.11.1 Alternative name(s): Cell death protein RIP Receptor-interacting protein 1 Short name= RIP-1 Serine/threonine-protein kinase RIP
Expression Region	1-656aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	78.9 kDa
Protein Length	Full Length



Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.