





Recombinant Vaccinia virus Protein OPG154 (OPG154)

Product Code	CSB-BP324948VAA
Abbreviation	Recombinant Vaccinia virus OPG154 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P20535
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Vaccinia virus (strain Copenhagen) (VACV)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	MDGTLFPGDDDLAIPATEFFSTKADKKPEAKREAIVKADEDDNEETLKQRLTNL EKKITNVTTKFEQIEKCKRNDEVLFRLENHAETLRAAMISLAKKIDVQTGRRPY E
Research Area	Others
Source	Baculovirus
Target Names	OPG154
Protein Names	Recommended name: 14 kDa fusion protein
Expression Region	1-110aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	16.5
Protein Length	Full Length

Image




Based on the SEQUEST from database of Baculovirus host and target protein, the LC-MS/MS Analysis result of CSB-BP324948VAA could indicate that this peptide derived from Baculovirus-expressed Vaccinia virus (strain Copenhagen) (VACV) A27L.



Based on the SEQUEST from database of Baculovirus host and target protein, the LC-MS/MS Analysis result of CSB-BP324948VAA could indicate that this peptide derived from Baculovirus-expressed Vaccinia virus (strain



Copenhagen) (VACV) A27L.



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.