



Myc Antibody, HRP conjugated

Product Code	CSB-PA015270EB01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01108
Immunogen	Recombinant Mouse Myc proto-oncogene protein (1-439AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Myc proto-oncogene protein (Proto-oncogene c-Myc) (Transcription factor p64), Myc
Immunogen Species	Mus musculus (Mouse)
Research Area	Epigenetics and Nuclear Signaling
Target Names	Myc
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.