



CAMK2N1 Antibody, Biotin conjugated

Product Code	CSB-PA004469ND01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7Z7J9
Immunogen	Peptide sequence from Human Calcium/calmodulin-dependent protein kinase II inhibitor 1 protein (9-27AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Calcium/calmodulin-dependent protein kinase II inhibitor 1 (CaMKII inhibitory protein alpha) (CaMKIIN-alpha), CAMK2N1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	CAMK2N1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.