



CHML Antibody, HRP conjugated

Product Code	CSB-PA005357LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P26374
Immunogen	Recombinant Human Rab proteins geranylgeranyltransferase component A 2 protein (480-656AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Rab proteins geranylgeranyltransferase component A 2 (Choroideremia-like protein) (Rab escort protein 2) (REP-2), CHML, REP2
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	CHML
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.