



EPHA8 Antibody, FITC conjugated

Product Code	CSB-PA007728LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P29322
Immunogen	Recombinant Human Ephrin type-A receptor 8 protein (45-249AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ephrin type-A receptor 8 (EC 2.7.10.1) (EPH- and ELK-related kinase) (EPH-like kinase 3) (EK3) (hEK3) (Tyrosine-protein kinase receptor EEK), EPHA8, EEK HEK3 KIAA1459
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	EPHA8
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.