



EXOC8 Antibody, FITC conjugated

Product Code	CSB-PA007887LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IYI6
Immunogen	Recombinant Human Exocyst complex component 8 protein (508-725AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Exocyst complex component 8 (Exocyst complex 84 kDa subunit), EXOC8
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	EXOC8
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.