



# EXOC8 Antibody, Biotin conjugated

Product Code	CSB-PA007887LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IYI6
Immunogen	Recombinant Human Exocyst complex component 8 protein (508-725AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Exocyst complex component 8 (Exocyst complex 84 kDa subunit), EXOC8
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	EXOC8