



HSPD1 Antibody, Biotin conjugated

Product Code	CSB-PA010847LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P10809
Immunogen	Recombinant Human 60 kDa heat shock protein, mitochondrial protein (225-573AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	60 kDa heat shock protein, mitochondrial (EC 3.6.4.9) (60 kDa chaperonin) (Chaperonin 60) (CPN60) (Heat shock protein 60) (HSP-60) (Hsp60) (HuCHA60) (Mitochondrial matrix protein P1) (P60 lymphocyte protein), HSPD1, HSP60
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	HSPD1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.