



TMEM131L Antibody, Biotin conjugated

Product Code	CSB-PA012218ND01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A2VDJ0
Immunogen	Peptide sequence from Human Transmembrane protein 131-like protein (963-978AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transmembrane protein 131-like, TMEM131L, KIAA0922
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TMEM131L