



LGALS4 Antibody, HRP conjugated

Product Code	CSB-PA012889LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P56470
Immunogen	Recombinant Human Galectin-4 protein (1-323AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Galectin-4 (Gal-4) (Antigen NY-CO-27) (L-36 lactose-binding protein) (L36LBP) (Lactose-binding lectin 4), LGALS4
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	LGALS4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.