



LILRB5 Antibody, HRP conjugated

Product Code	CSB-PA012943LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75023
Immunogen	Recombinant Human Leukocyte immunoglobulin-like receptor subfamily B member 5 protein (24-350AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Leukocyte immunoglobulin-like receptor subfamily B member 5 (CD85 antigen-like family member C) (Leukocyte immunoglobulin-like receptor 8) (LIR-8) (CD antigen CD85c), LILRB5, LIR8
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	LILRB5