



MAP4K4 Antibody, HRP conjugated

Product Code	CSB-PA013439LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95819
Immunogen	Recombinant Human Mitogen-activated protein kinase kinase kinase kinase 4 protein (866-1104AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mitogen-activated protein kinase kinase kinase kinase 4 (EC 2.7.11.1) (HPK/GCK-like kinase HGK) (MAPK/ERK kinase kinase kinase 4) (MEK kinase kinase 4) (MEKKK 4) (Nck-interacting kinase), MAP4K4, HGK KIAA0687 NIK
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MAP4K4