



MAP4K4 Antibody, Biotin conjugated

Product Code	CSB-PA013439LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95819
Immunogen	Recombinant Human Mitogen-activated protein kinase kinase kinase kinase 4 protein (866-1104AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mitogen-activated protein kinase kinase kinase kinase 4 (EC 2.7.11.1) (HPK/GCK-like kinase HGK) (MAPK/ERK kinase kinase kinase 4) (MEK kinase kinase 4) (MEKKK 4) (Nck-interacting kinase), MAP4K4, HGK KIAA0687 NIK
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MAP4K4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.