



MFF Antibody, FITC conjugated

Product Code	CSB-PA013751LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9GZY8
Immunogen	Recombinant Human Mitochondrial fission factor protein (1-322AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mitochondrial fission factor, MFF, C2orf33
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	MFF
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.