



LPAR6 Antibody, HRP conjugated

Product Code	CSB-PA017336OB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P43657
Immunogen	Peptide sequence from Human Lysophosphatidic acid receptor 6 protein (292-313AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Lysophosphatidic acid receptor 6 (LPA receptor 6) (LPA-6) (Oleoyl-L-alpha-lysophosphatidic acid receptor) (P2Y purinoceptor 5) (P2Y5) (Purinergic receptor 5) (RB intron encoded G-protein coupled receptor), LPAR6, P2RY5
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	LPAR6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.