



PITPNA Antibody, Biotin conjugated

Product Code	CSB-PA018035LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q00169
Immunogen	Recombinant Human Phosphatidylinositol transfer protein alpha isoform protein (1-270AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Phosphatidylinositol transfer protein alpha isoform (PI-TP-alpha) (PtdIns transfer protein alpha) (PtdInsTP alpha), PITPNA, PITPN
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PITPNA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.