



PLG Antibody, Biotin conjugated

Product Code	CSB-PA018188LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P00747
Immunogen	Recombinant Human Plasminogen protein (103-343AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Plasminogen (EC 3.4.21.7) [Cleaved into: Plasmin heavy chain A; Activation peptide; Angiostatin; Plasmin heavy chain A, short form; Plasmin light chain B], PLG
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	PLG
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.