



PON1 Antibody, Biotin conjugated

Product Code	CSB-PA018369LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P27169
Immunogen	Recombinant Human Serum paraoxonase/arylesterase 1 protein (103-355AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serum paraoxonase/arylesterase 1 (PON 1) (EC 3.1.1.2) (EC 3.1.1.81) (EC 3.1.8.1) (Aromatic esterase 1) (A-esterase 1) (K-45) (Serum aryldialkylphosphatase 1), PON1, PON
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	PON1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.