



PRAMEF10 Antibody, FITC conjugated

Product Code	CSB-PA018605LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60809
Immunogen	Recombinant Human PRAME family member 10 protein (1-261AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	PRAME family member 10, PRAMEF10
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	PRAMEF10