



S100A7 Antibody, Biotin conjugated

Product Code	CSB-PA020635HD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P31151
Immunogen	Recombinant Human Protein S100-A7 protein (2-101AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein S100-A7 (Psoriasin) (S100 calcium-binding protein A7), S100A7, PSOR1 S100A7C
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	S100A7
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.