



ST5 Antibody, HRP conjugated

Product Code	CSB-PA022758NB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P78524
Immunogen	Peptide sequence from Human Suppression of tumorigenicity 5 protein (68-88AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Suppression of tumorigenicity 5 protein (DENN domain-containing protein 2B) (HeLa tumor suppression 1), ST5, DENND2B HTS1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ST5