



ST5 Antibody, Biotin conjugated

Product Code	CSB-PA022758ND01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P78524
Immunogen	Peptide sequence from Human Suppression of tumorigenicity 5 protein (68-88AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Suppression of tumorigenicity 5 protein (DENN domain-containing protein 2B) (HeLa tumor suppression 1), ST5, DENND2B HTS1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ST5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.