



TJP1 Antibody, HRP conjugated

| | |
|----------------------------|---|
| Product Code | CSB-PA023573OB01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q07157 |
| Immunogen | Peptide sequence from Human Tight junction protein ZO-1 protein (116-138AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | HRP |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | Antigen Affinity Purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Tight junction protein ZO-1 (Tight junction protein 1) (Zona occludens protein 1) (Zonula occludens protein 1), TJP1, ZO1 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Tags & Cell Markers |
| Target Names | TJP1 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |