



TJP1 Antibody, FITC conjugated

Product Code	CSB-PA023573OC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q07157
Immunogen	Peptide sequence from Human Tight junction protein ZO-1 protein (116-138AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tight junction protein ZO-1 (Tight junction protein 1) (Zona occludens protein 1) (Zonula occludens protein 1), TJP1, ZO1
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	TJP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.