



EMC4 Antibody, Biotin conjugated

Product Code	CSB-PA023887OD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5J8M3
Immunogen	Peptide sequence from Human ER membrane protein complex subunit 4 protein (43-60AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	ER membrane protein complex subunit 4 (Cell proliferation-inducing gene 17 protein) (Transmembrane protein 85), EMC4, TMEM85
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	EMC4