



TNFRSF1A Antibody, Biotin conjugated

Product Code	CSB-PA023977LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P19438
Immunogen	Recombinant Human Tumor necrosis factor receptor superfamily member 1A protein (249-455AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tumor necrosis factor receptor superfamily member 1A (Tumor necrosis factor receptor 1) (TNF-R1) (Tumor necrosis factor receptor type I) (TNF-RI) (TNFR-I) (p55) (p60) (CD antigen CD120a) [Cleaved into: Tumor necrosis factor receptor superfamily member 1A, membrane form; Tumor necrosis factor-binding protein 1 (TBPI)], TNFRSF1A, TNFAR TNFR1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TNFRSF1A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.