



CD24 Antibody, Biotin conjugated

Product Code	CSB-PA10468D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P25063
Immunogen	Peptide sequence from Human Signal transducer CD24 protein (28-46AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Signal transducer CD24 (Small cell lung carcinoma cluster 4 antigen) (CD antigen CD24), CD24, CD24A
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	CD24
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.