




Uncharacterized protein Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA27779A0Rb |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | E7EVH7 |
| Immunogen | Recombinant Human Uncharacterized protein (502-717AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA, WB; Recommended dilution: WB:1:500-1:2000 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | Antigen Affinity Purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Uncharacterized protein |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Others |
| Target Names | Uncharacterized protein |

Image



Western Blot

Positive WB detected in: HepG2 whole cell lysate, 293T whole cell lysate, SH-SY5Y whole cell lysate

All lanes: Uncharacterized protein antibody at 1:1000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 84 kDa

Observed band size: 84 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.