




Uncharacterized protein Antibody

Product Code	CSB-PA27779A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	E7EVH7
Immunogen	Recombinant Human Uncharacterized protein (502-717AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Uncharacterized protein
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	Uncharacterized protein

Image



Western Blot

Positive WB detected in: HepG2 whole cell lysate, 293T whole cell lysate, SH-SY5Y whole cell lysate

All lanes: Uncharacterized protein antibody at 1:1000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 84 kDa

Observed band size: 84 kDa