



TRAPPC8 Antibody, Biotin conjugated

Product Code	CSB-PA27999D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y2L5
Immunogen	Recombinant Human Trafficking protein particle complex subunit 8 protein (1205-1435AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Trafficking protein particle complex subunit 8 (Protein TRS85 homolog), TRAPPC8, KIAA1012
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TRAPPC8
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.