



VP2 Antibody, Biotin conjugated

Product Code	CSB-PA302471LD01CID
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P69484
Immunogen	Recombinant Chicken anemia virus Dual specificity protein phosphatase VP2 protein (1-216AA)
Raised In	Rabbit
Species Reactivity	Chicken anemia virus
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Dual specificity protein phosphatase VP2 (EC 3.1.3.16) (EC 3.1.3.48), VP2
Immunogen Species	Chicken anemia virus (isolate Germany Cuxhaven-1) (CAV)
Research Area	Others
Target Names	VP2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.