




N Antibody

Product Code	CSB-PA339612HA01BKT
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P35943
Immunogen	Recombinant Bovine respiratory syncytial virus Nucleoprotein protein (1-391AA)
Raised In	Rabbit
Species Reactivity	Bovine respiratory syncytial virus
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nucleoprotein (Protein N) (Nucleocapsid protein), N
Immunogen Species	Bovine respiratory syncytial virus (strain 391-2) (BRS)
Research Area	Others
Target Names	N

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: N antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 59 kDa

Observed band size: 59 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.