



N Antibody, Biotin conjugated

Product Code	CSB-PA339612HD01BKT
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P35943
Immunogen	Recombinant Bovine respiratory syncytial virus Nucleoprotein protein (1-391AA)
Raised In	Rabbit
Species Reactivity	Bovine respiratory syncytial virus
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nucleoprotein (Protein N) (Nucleocapsid protein), N
Immunogen Species	Bovine respiratory syncytial virus (strain 391-2) (BRS)
Research Area	Others
Target Names	N