




M Antibody

Product Code	CSB-PA356390LA01VBG
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04876
Immunogen	Recombinant Vesicular stomatitis Indiana virus Matrix protein (1-237AA)
Raised In	Rabbit
Species Reactivity	Vesicular stomatitis Indiana virus
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Matrix protein, M
Immunogen Species	Vesicular stomatitis Indiana virus (strain Glasgow) (VSIV)
Research Area	Others
Target Names	M

Image



Western Blot

Positive WB detected in Recombinant protein
All lanes: M antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution
Predicted band size: 35 kDa
Observed band size: 35 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.