



M Antibody, Biotin conjugated

Product Code	CSB-PA356390LD01VBG
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04876
Immunogen	Recombinant Vesicular stomatitis Indiana virus Matrix protein (1-237AA)
Raised In	Rabbit
Species Reactivity	Vesicular stomatitis Indiana virus
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Matrix protein, M
Immunogen Species	Vesicular stomatitis Indiana virus (strain Glasgow) (VSIV)
Research Area	Others
Target Names	M
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.