



HSPA1A Antibody, Biotin conjugated

Product Code	CSB-PA357302LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0DMV8
Immunogen	Recombinant Human Heat shock 70 kDa protein 1A protein (291-641AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Heat shock 70 kDa protein 1A (Heat shock 70 kDa protein 1) (HSP70-1) (HSP70.1), HSPA1A, HSP72 HSPA1 HSX70
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	HSPA1A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.