



SHISA9 Antibody, HRP conjugated

Product Code	CSB-PA457545LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	B4DS77
Immunogen	Recombinant Human Protein shisa-9 protein (171-424AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein shisa-9, SHISA9
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	SHISA9
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.