




algL Antibody

Product Code	CSB-PA524422LA01DPX
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O52195
Immunogen	Recombinant Azotobacter vinelandii Alginate lyase protein (24-374AA)
Raised In	Rabbit
Species Reactivity	Azotobacter vinelandii
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Alginate lyase (EC 4.2.2.3) (Poly(beta-D-mannuronate) lyase), algL
Immunogen Species	Azotobacter vinelandii
Research Area	Others
Target Names	algL

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: algL antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 62 kDa

Observed band size: 62 kDa