



# PAP1 Antibody, FITC conjugated

<b>Product Code</b>	CSB-PA529312LC01DOA
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	O81439
<b>Immunogen</b>	Recombinant Arabidopsis thaliana Probable plastid-lipid-associated protein 1, chloroplastic protein (56-318AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Arabidopsis thaliana
<b>Form</b>	Liquid
<b>Conjugate</b>	FITC
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Probable plastid-lipid-associated protein 1, chloroplastic (AtPap1) (Fibrillin-1a) (Plastoglobulin 35) (AtPGL35), PAP1, FBN1a FIB1a PGL35
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Research Area</b>	Cardiovascular
<b>Target Names</b>	PAP1
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.