



PAP1 Antibody, Biotin conjugated

Product Code	CSB-PA529312LD01DOA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O81439
Immunogen	Recombinant Arabidopsis thaliana Probable plastid-lipid-associated protein 1, chloroplastic protein (56-318AA)
Raised In	Rabbit
Species Reactivity	Arabidopsis thaliana
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Probable plastid-lipid-associated protein 1, chloroplastic (AtPap1) (Fibrillin-1a) (Plastoglobulin 35) (AtPGL35), PAP1, FBN1a FIB1a PGL35
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Research Area	Cardiovascular
Target Names	PAP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.