



FABP9 Antibody, Biotin conjugated

Product Code	CSB-PA604514HD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q0Z7S8
Immunogen	Recombinant Human Fatty acid-binding protein 9 protein (1-132AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Fatty acid-binding protein 9 (Testis lipid-binding protein) (TLBP) (Testis-type fatty acid-binding protein) (T-FABP), FABP9
Immunogen Species	Homo sapiens (Human)
Research Area	Developmental Biology
Target Names	FABP9
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.