



RAB39A Antibody, Biotin conjugated

Product Code	CSB-PA614818LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14964
Immunogen	Recombinant Human Ras-related protein Rab-39A protein (1-217AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ras-related protein Rab-39A (Rab-39), RAB39A, RAB39
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	RAB39A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.